

## PRPF4 RABBIT MAB

**Cat.#:** N262771

**Product Name:** Anti-PRPF4 Rabbit Monoclonal Antibody

**Synonyms:** PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60

**UNIPROT ID:** O43172

**Background:** Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

**Immunogen:** A synthetic peptide of human PRPF4

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R07-6C5

**MW:** Calculated MW: 58 kDa; Observed MW: 58 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---