

PRPF4 RABBIT MAB

Cat.#: N262771

Product Name: Anti-PRPF4 Rabbit Monoclonal Antibody

Synonyms: PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60

UNIPROT ID: O43172

Background: Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

Immunogen: A synthetic peptide of human PRPF4

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-6C5

MW: Calculated MW: 58 kDa; Observed MW: 58 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
