

RAB25 RABBIT PAB

Cat.#: S219626

Product Name: Anti-RAB25 Rabbit Polyclonal Antibody

Synonyms: CATX-8, RAB11C

UNIPROT ID: P57735 (Gene Accession - NP_065120)

Background: RAB proteins, such as RAB25, are members of the RAS superfamily of small GTPases that are involved in membrane trafficking. Members of the RAB11 subfamily, including RAB25, control the return of internalized membrane-associated moieties to the cell surface.

Immunogen: Synthetic peptide of human RAB25

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-150; ELISA: 1000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Signal Transduction, Cancer

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver cancer tissue using 219626(RAB25 Antibody) at a dilution of 1/50(Cell membrane, Cytoplasm).