

RB RABBIT MAB

Cat.#: N261697

Product Name: Anti-Rb Rabbit Monoclonal Antibody

Synonyms: p110-RB1; pRb; pp105

UNIPROT ID: P13405

Background: Tumor suppressor that is a key regulator of the G1/S transition of the cell cycle.

Immunogen: Recombinant protein of mouse Rb

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-8G7

MW: Calculated MW: 105 kDa; Observed MW: 105 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Rb in C6, Jurkat, HeLa lysates using Rb antibody.