

## RCOR3 RABBIT PAB

**Cat.#:** S218535

**Product Name:** Anti-RCOR3 Rabbit Polyclonal Antibody

**Synonyms:**

**UNIPROT ID:** Q9P2K3 (Gene Accession - BC031608 )

**Background:** May act as a component of a corepressor complex that represses transcription.

**Immunogen:** Fusion protein of human RCOR3

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 20-100; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 218535(RCOR3 Antibody) at a dilution of 1/20(Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the fusion protein and then with 218535(Anti-RCOR3 Antibody) at dilution 1/20.