

RPA70 RABBIT MAB

Cat.#: N262846

Product Name: Anti-RPA70 Rabbit Monoclonal Antibody

Synonyms: RPA1; HSSB; MST075; REPA1; RF-A; RP-A; RPA70

UNIPROT ID: P27694

Background: Plays an essential role in several cellular processes in DNA metabolism including replication, recombination and DNA repair. Binds and subsequently stabilizes single-stranded DNA intermediates and thus prevents complementary DNA from reannealing.

Immunogen: A synthetic peptide of human RPA70

Applications: WB,IHC-F,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-7B7

MW: Calculated MW: 68 kDa; Observed MW: 68 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RPA70 in A549 lysates using RPA70 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
