

## RRM1 RABBIT PAB

**Cat.#:** N225989

**Product Name:** Anti-RRM1 Rabbit pAb

**Synonyms:** RRM1; R1; RIR1; RR1

**UNIPROT ID:** P23921

**Background:** Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the corresponding ribonucleotides.

**Immunogen:** A synthesized peptide derived from human RRM1

**Applications:** WB,IHC-P,ICC/IF,IP,FC

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 90 kDa; Observed MW: 90 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RRM1 in A549 lysates using RRM1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---