

S100 ALPHA9 RABBIT MAB

Cat.#: N262860

Product Name: Anti-S100 alpha9 Rabbit Monoclonal Antibody

Synonyms: S100 alpha9; Calgranulin-B; Leukocyte L1 complex heavy chain

UNIPROT ID: P31725

Background: S100A9 is a calcium- and zinc-binding protein which plays a prominent role in the regulation of inflammatory processes and immune response.

Immunogen: Recombinant protein of mouse S100A9

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-3C3

MW: Calculated MW: 13 kDa; Observed MW: 13 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of S100A9 in mouse spleen, mouse lung lysates using S100 alpha9 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
