

## SCGB1A1 RABBIT PAB

**Cat.#:** S216873

**Product Name:** Anti-SCGB1A1 Rabbit Polyclonal Antibody

**Synonyms:** UGB; UPI; CC10; CC16; CCSP; UP-1; CCPBP

**UNIPROT ID:** P11684 (Gene Accession - BC004481 )

**Background:** This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.

**Immunogen:** Fusion protein of human SCGB1A1

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 50-200; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using 216873(SCGB1A1 Antibody) at a dilution of 1/80(Cytoplasm).