



SLC15A1 RABBIT PAB**Cat.#:** S221775**Product Name:** Anti-SLC15A1 Rabbit Polyclonal Antibody**Synonyms:** PEPT1; HPECT1; HPEPT1**UNIPROT ID:** P46059 (Gene Accession - NP_005064)**Background:** This gene encodes an intestinal hydrogen peptide cotransporter that is a member of the solute carrier family 15. The encoded protein is localized to the brush border membrane of the intestinal epithelium and mediates the uptake of di- and tripeptides from the lumen into the enterocytes. This protein plays an important role in the uptake and digestion of dietary proteins. This protein also facilitates the absorption of numerous peptidomimetic drugs.**Immunogen:** Synthetic peptide of human SLC15A1**Applications:** ELISA, IHC**Recommended Dilutions:** IHC: Oct-50; ELISA: 5000-10000**Host Species:** Rabbit**Clonality:** Rabbit Polyclonal**Isotype:** Immunogen-specific rabbit IgG**Purification:** Antigen affinity purification**Species Reactivity:** Human**Constituents:** PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol**Research Areas:** Metabolism**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing

Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 221775(SLC15A1 Antibody) at a dilution of 1/20(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the synthetic peptide and then with 221775(Anti-SLC15A1 Antibody) at dilution 1/20.