

SNX4 RABBIT MAB

Cat.#: N262922

Product Name: Anti-SNX4 Rabbit Monoclonal Antibody

Synonyms: ATG24B

UNIPROT ID: O95219

Background: Involved in the regulation of endocytosis and in several stages of intracellular trafficking (PubMed:12668730). Plays a role in recycling endocytosed transferrin receptor and prevent its degradation.

Immunogen: A synthetic peptide of human SNX4

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-8A9

MW: Calculated MW: 52 kDa; Observed MW: 52 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SNX4 in C6, 3T3, Hela lysates using SNX4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
