

SOX17 RABBIT PAB

Cat.#: N225160

Product Name: Anti-SOX17 Rabbit pAb

Synonyms: SOX17; Transcription factor SOX-17

UNIPROT ID: Q9H6I2

Background: Acts as transcription regulator that binds target promoter DNA and bends the DNA.

Immunogen: Synthesized peptide derived from Sox-17 . at AA range: 90-170

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 44 kDa; Observed MW: 44 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SOX17 in Jurkat lysates using Sox17 antibody.