

## STAR RABBIT PAB

**Cat.#:** N225884

**Product Name:** Anti-StAR Rabbit pAb

**Synonyms:** STARD1

**UNIPROT ID:** P49675

**Background:** Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer membrane to the inner mitochondrial membrane.

**Immunogen:** A synthesized peptide derived from human StAR

**Applications:** WB,IHC-P

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 32 kDa; Observed MW: 32 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of StAR in SKBR3 lysates using StAR antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---