

STAT6 RABBIT MAB

Cat.#: N263563

Product Name: Anti-STAT6 Rabbit Monoclonal Antibody

Synonyms: STAT6; Signal transducer and activator of transcription 6; IL-4 Stat

UNIPROT ID: P42226

Background: STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

Immunogen: A synthetic peptide of human STAT6

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-0B5

MW: Calculated MW: 94 kDa; Observed MW: 110 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STAT6 in HeLa, U2OS lysates using STAT6 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
