

## STRIATIN 4 RABBIT MAB

**Cat.#:** N262948

**Product Name:** Anti-Striatin 4 Rabbit Monoclonal Antibody

**Synonyms:** ZIN; PPP2R6C; zinedin

**UNIPROT ID:** Q9NRL3

**Background:** Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein. Miscellaneous The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).

**Immunogen:** Recombinant protein of human Striatin 4

**Applications:** WB, ICC/IF, IP

**Recommended Dilutions:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R04-4J1

**MW:** Calculated MW: 81 kDa; Observed MW: 81 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Striatin 4 in Hela lysates using Striatin 4 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---