

STRIATIN 4 RABBIT MAB

Cat.#: N262948

Product Name: Anti-Striatin 4 Rabbit Monoclonal Antibody

Synonyms: ZIN; PPP2R6C; zinedin

UNIPROT ID: Q9NRL3

Background: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein. Miscellaneous The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).

Immunogen: Recombinant protein of human Striatin 4

Applications: WB, ICC/IF, IP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-4J1

MW: Calculated MW: 81 kDa; Observed MW: 81 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Striatin 4 in Hela lysates using Striatin 4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
