

SYNDECAN 3 RABBIT MAB

Cat.#: N262959

Product Name: Anti-Syndecan 3 Rabbit Monoclonal Antibody

Synonyms: SDCN; SYND3

UNIPROT ID: O75056

Background: Cell surface proteoglycan that may bear heparan sulfate . May have a role in the organization of cell shape by affecting the actin cytoskeleton, possibly by transferring signals from the cell surface in a sugar-dependent mechanism.

Immunogen: A synthetic peptide of human Syndecan 3

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-6B2

MW: Calculated MW: 45 kDa; Observed MW: 60 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Syndecan 3 in K562, C6, 3T3 lysates using Syndecan 3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
