

## T BOX PROTEIN 3 RABBIT MAB

**Cat.#:** N262976

**Product Name:** Anti-T Box Protein 3 Rabbit Monoclonal Antibody

**Synonyms:** T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL

**UNIPROT ID:** O15119

**Background:** Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

**Immunogen:** A synthetic peptide of human Tbx3

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R07-5B9

**MW:** Calculated MW: 79 kDa; Observed MW: 79 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Mouse


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.