

T BOX PROTEIN 3 RABBIT MAB

Cat.#: N262976

Product Name: Anti-T Box Protein 3 Rabbit Monoclonal Antibody

Synonyms: T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL

UNIPROT ID: O15119

Background: Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

Immunogen: A synthetic peptide of human Tbx3

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-5B9

MW: Calculated MW: 79 kDa; Observed MW: 79 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.