

TAGAP RABBIT MAB

Cat.#: N262967

Product Name: Anti-TAGAP Rabbit Monoclonal Antibody

Synonyms: FKSG15; IDDM21; TAGAP1; ARHGAP47

UNIPROT ID: Q8N103

Background: May function as a GTPase-activating protein and may play important roles during T-cell activation.

Immunogen: Recombinant protein of human TAGAP

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-3D3

MW: Calculated MW: 81 kDa; Observed MW: 81 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TAGAP in rat Brain, C6, CHO-K1, HeLa lysates using TAGAP antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
