

TCEAL6 RABBIT PAB

Cat.#: N225200

Product Name: Anti-TCEAL6 Rabbit pAb

Synonyms: TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6

UNIPROT ID: Q6IPX3

Background: May be involved in transcriptional regulation.

Immunogen: Synthesized peptide derived from TCEAL6 . at AA range: 1-80

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 22 kDa; Observed MW: 32 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TCEAL6 in CoLo lysates using TCEAL6 antibody.