

## TCTP RABBIT MAB

**Cat.#:** N262981

**Product Name:** Anti-TCTP Rabbit Monoclonal Antibody

**Synonyms:** Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23

**UNIPROT ID:** P13693

**Background:** Involved in calcium binding and microtubule stabilization.

**Immunogen:** A synthetic peptide of human TCTP

**Applications:** WB, ICC/IF, IP

**Recommended Dilutions:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R02-5D7

**MW:** Calculated MW: 20 kDa; Observed MW: 20 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human, Mouse, Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of TCTP (green) in HeLa using TCTP antibody, and DAPI (blue)



Western blot analysis of TCTP in K562, rat Brain, C6, 3T3, HeLa lysates using TCTP antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---