

## TCTP RABBIT PAB

**Cat.#:** N225857

**Product Name:** Anti-TCTP Rabbit pAb

**Synonyms:** Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23

**UNIPROT ID:** P13693

**Background:** Involved in calcium binding and microtubule stabilization.

**Immunogen:** A synthesized peptide derived from human TCTP

**Applications:** WB,IHC-P,ICC/IF,IP,FC

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 20 kDa; Observed MW: 20 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---