

TMEFF1 RABBIT PAB

Cat.#: S217930

Product Name: Anti-TMEFF1 Rabbit Polyclonal Antibody

Synonyms: TR-1; H7365; C9orf2; CT120.1

UNIPROT ID: Q8IYR6 (Gene Accession - BC035056)

Background: May inhibit NODAL and BMP signaling during neural patterning (By similarity). May be a tumor suppressor in brain cancers.

Immunogen: Fusion protein of human TMEFF1

Applications: ELISA, IHC

Recommended Dilutions: IHC: 150-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Cancer

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 217930(TMEFF1 Antibody) at a dilution of 1/140(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 217930(Anti-TMEFF1 Antibody) at dilution 1/140.