

TNFAIP8 RABBIT MAB

Cat.#: N261564

Product Name: Anti-TNFAIP8 Rabbit Monoclonal Antibody

Synonyms: TNF alpha-induced protein 8; NDED; SCC-S2

UNIPROT ID: O95379

Background: Acts as a negative mediator of apoptosis and may play a role in tumor progression.

Immunogen: A synthetic peptide of human TNFAIP8

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-5F3

MW: Calculated MW: 23 kDa; Observed MW: 21-23 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TNFAIP8 in CEM, K562 lysates using TNFAIP8 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
