

TNFAIP8L3 RABBIT PAB

Cat.#: N225228

Product Name: Anti-TNFAIP8L3 Rabbit pAb

Synonyms: Tumor necrosis factor alpha-induced protein 8-like protein 3; TNF alpha-induced protein 8-like protein 3; TNFAIP8-like protein 3

UNIPROT ID: Q5GJ75

Background: Acts as a lipid transfer protein.

Immunogen: Synthesized peptide derived from human TNFAIP8L3. at AA range: 15-45

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 33 kDa; Observed MW: 33 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TNFAIP8L3 in 293T , mouse brain lysates using TNFAIP8L3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
