

## TRANSMEMBRANE PROTEIN 72 RABBIT PAB

**Cat.#:** N225244

**Product Name:** Anti-Transmembrane Protein 72 Rabbit pAb

**Synonyms:** Transmembrane protein 72; Kidney-specific secretory protein of 37 kDa; C10orf127; KSP37

**UNIPROT ID:** A0PK05

**Background:** This gene encodes a transmembrane protein which may be expressed specifically in the kidney.

**Immunogen:** Synthesized peptide derived from human Tmem72  
Polyclonal

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 30 kDa; Observed MW: 30 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Transmembrane Protein 72 in various lysates using Transmembrane Protein 72 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---