

TRIM47 RABBIT PAB

Cat.#: S214827

Product Name: Anti-TRIM47 Rabbit Polyclonal Antibody

Synonyms: GOA; RNF100

UNIPROT ID: Q96LD4 (Gene Accession - NP_258411)

Background: Enables ubiquitin-protein transferase activity. Involved in protein ubiquitination. Located in cytosol. [provided by Alliance of Genome Resources, Apr 2022]

Immunogen: Synthetic peptide of human TRIM47

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-200; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cancer

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 214827 (TRIM47 Antibody) at a dilution of 1/55 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the synthetic peptide and then with 214827 (Anti-TRIM47 Antibody) at dilution 1/55.