

## TRYPTASE DELTA RABBIT PAB

**Cat.#:** N225019

**Product Name:** Anti-Tryptase delta Rabbit pAb

**Synonyms:** Tryptase delta (EC 3.4.21.59; Delta-tryptase; HmMCP-3-like tryptase III; Mast cell mMCP-7-like; Tryptase-3)

**UNIPROT ID:** Q9BZJ3

**Background:** Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type.

**Immunogen:** Synthesized peptide derived from Tryptase-3 at AA range: 51-100

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 27 kDa; Observed MW: 27 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** -

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tryptase delta in HEPG2 lysates using Tryptase delta antibody.