

TRYPTASE DELTA RABBIT PAB

Cat.#: N225019

Product Name: Anti-Tryptase delta Rabbit pAb

Synonyms: Tryptase delta (EC 3.4.21.59; Delta-tryptase; HmMCP-3-like tryptase III; Mast cell mMCP-7-like; Tryptase-3)

UNIPROT ID: Q9BZJ3

Background: Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type.

Immunogen: Synthesized peptide derived from Tryptase-3 at AA range: 51-100

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 27 kDa; Observed MW: 27 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tryptase delta in HEPG2 lysates using Tryptase delta antibody.