

## USP38 RABBIT PAB

**Cat.#:** S218330

**Product Name:** Anti-USP38 Rabbit Polyclonal Antibody

**Synonyms:** HP43.8KD

**UNIPROT ID:** Q8NB14 (Gene Accession - BC068975 )

**Background:** USP38 (ubiquitin specific peptidase 38), also known as HP43.8KD, is a 1,042 amino acid deubiquitinating enzyme that participates in the Ub pathway, exhibiting a preference towards 'Lys-63'-linked Ubiquitin chains. A member of the peptidase C19 family, the catalytic activity of USP38 involves a combination of the ubiquitin carboxyl-terminal thiolester and water to produce ubiquitin and a thiol. USP38 is highly expressed in skeletal muscle and adrenal gland.

**Immunogen:** Fusion protein of human USP38

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 25-100; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human, Mouse


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Epigenetics and Nuclear Signaling, Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human tonsil tissue using 218330(USP38 Antibody) at a dilution of 1/30(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with the fusion protein and then with 218330(Anti-USP38 Antibody) at dilution 1/30.