

## VEGF RECEPTOR 1 RABBIT PAB

**Cat.#:** N225566

**Product Name:** Anti-VEGF Receptor 1 Rabbit pAb

**Synonyms:** VEGFR-1; VEGF Receptor 1; FLT-1; FLT; FRT; VEGFR1; Fms-like tyrosine kinase 1; Tyrosine-protein kinase FRT

**UNIPROT ID:** P17948

**Background:** Tyrosine-protein kinase that acts as a cell-surface receptor for VEGFA, VEGFB and PGF, and plays an essential role in the development of embryonic vasculature, the regulation of angiogenesis, cell survival, cell migration, macrophage function, chemotaxis, and cancer cell invasion. May play an essential role as a negative regulator of embryonic angiogenesis by inhibiting excessive proliferation of endothelial cells.

**Immunogen:** A synthesized peptide derived from human VEGF Receptor 1

**Applications:** WB,IHC-P,ICC/IF,IP,FC

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 151 kDa; Observed MW: 151 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cardiovascular

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of VEGFR1 in A431 lysates using VEGF Receptor 1 antibody.