

VPS33A RABBIT PAB

Cat.#: S218015

Product Name: Anti-VPS33A Rabbit Polyclonal Antibody

Synonyms: MPSPS

UNIPROT ID: Q96AX1 (Gene Accession - BC016617)

Background: This gene encodes a tethering protein and a core subunit of the homotypic fusion and protein sorting (HOPS) complex. The HOPS complex and a second endosomal tethering complex called the class C core vacuole/endosome tethering (CORVET) complex perform diverse functions in endocytosis including membrane tethering, RabGTPase interaction, activation and proofreading of synaptic-soluble N-ethylmaleimide-sensitive factor attachment receptor (SNARE) assembly to drive membrane fusion, and endosome-to-cytoskeleton attachment. The HOPS complex controls endosome maturation as well as endosome traffic to the lysosome. This complex is essential for vacuolar fusion and is required for adaptor protein complex 3-dependent transport from the golgi to the vacuole. The encoded protein belongs to the Sec1/Munc18 (SM) family of SNARE-mediated membrane fusion regulators. Naturally occurring mutations in this gene are associated with a novel mucopolysaccharidosis-like disease.

Immunogen: Fusion protein of human VPS33A

Applications: ELISA, IHC

Recommended Dilutions: IHC: 100-200; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer tissue using 218015(VPS33A Antibody) at a dilution of 1/80(Cytoplasm and Cell membrane).