

ZHX2 RABBIT MAB

Cat.#: N261629

Product Name: Anti-ZHX2 Rabbit Monoclonal Antibody

Synonyms: AFP regulator 1; Zinc finger and homeodomain protein 2

UNIPROT ID: Q9Y6X8

Background: Acts as a transcriptional repressor.

Immunogen: Recombinant protein of human ZHX2

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-9B5

MW: Calculated MW: 92 kDa; Observed MW: 110 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of ZHX2 (green) in MCF-7 using ZHX2 antibody, and DAPI (blue)



Western blot analysis of ZHX2 in rat Brain, C6, HeLa lysates using ZHX2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
