

ZNF499 RABBIT PAB

Cat.#: N225519

Product Name: Anti-ZNF499 Rabbit pAb

Synonyms: Zinc finger protein 499; Zinc finger and BTB domain-containing protein 45; ZBTB45

UNIPROT ID: Q96K62

Background: May be involved in transcriptional regulation (Probable).

Immunogen: Recombinant Protein of ZBTB45

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 54 kDa; Observed MW: 70 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded mouse Kidney Tissue using ZNF499 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of ZNF499 in Hela, Jurkat, mouse Brain, mouse Kidney, rat Brain lysates using ZNF499 antibody.