


HA TAG MAB**Anti-HA Tag Mouse Monoclonal Antibody****Cat. #:** 28003**Size:** 100 µL**Symbol:** HA-tag**Description:** Anti-HA Tag Mouse Monoclonal Antibody**Background:** The HA-tag antibody recognizes any protein containing HA Tag.**Immunogen:** HA**Applications:** ELISA, WB, IF, IHC**Recommended Dilutions:**

ELISA 1:2000–1:5000

WB 1:1000–1:5000

IF 1:100–1:200

Concentration: 1 mg/ml**Host Species:** Mouse**Format:** Liquid**Clonality:** Monoclonal**Isotype:** IgG**Purity:** Purified from ascites**Constituents:** PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol**Species Reactivity:** Recognizes protein contain HA.**Storage Conditions:** Store at $-20^{\circ}C$. Avoid repeated freezing and thawing**Western blot:**

Western blot: Recombinant HA tag fusion protein was diluted with HA mouse monoclonal antibody (1:1000), and the protein loading amount was 50ng, and Western blot analysis was performed.

