

HA TAG MAB

Anti-HA Tag Mouse Monoclonal Antibody

Cat. #: 28003

Size: 100 μ L

Symbol: HA-tag

Description: Anti-HA Tag Mouse Monoclonal Antibody

Background: The HA-tag antibody recognizes any protein containing HA Tag.

Immunogen: HA

Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA 1:2000-1:5000

WB 1:1000-1:5000

IF 1:100-1:200

Concentration: 1 mg/ml

Host Species: Mouse

Format: Liquid

Clonality: Monoclonal

Isotype: IgG


Purity: Purified from ascites

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes protein contain HA.

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



Western blot: Recombinant HA tag fusion protein was diluted with HA mouse monoclonal antibody (1;1000), and the protein loading amount was 50ng, and Western blot analysis was performed.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
