


**APOLIPOPROTEIN A1****Apolipoprotein A1 Monoclonal Antibody****Cat. #:** 26025**Gene Symbol:** APOA1**Description:** Anti-Apolipoprotein A1 Mouse Monoclonal Antibody**Background:** Apolipoprotein A1 is the major protein component of high-density lipoprotein (HDL) ("good cholesterol") in plasma. Apolipoprotein A1 helps to clear cholesterol from arteries.**Immunogen:** Recombinant full length protein of human Apolipoprotein A1**Applications:** ELISA, WB, IHC, IP**Recommended Dilutions:**

ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:50–1:100

**Concentration:** 1 mg/ml**Host Species:** Mouse**Format:** Liquid**Clonality:** Monoclonal**Isotype:** IgG**Purity:** Purified from ascites**Preservative:** No**Constituents:** PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol**Species Reactivity:** Apolipoprotein A1 antibody recognizes Apolipoprotein A1 of primates.**Storage Conditions:** Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing**Western blot:**

WB: anti-Apo A1 mAb

**Western blot analysis of Apolipoprotein A1 proteins.**

