

BECKMAN COULTER BIOMEK I-SERIES 90UL PIPETTE TIPS

For use with Multichannel-96 Heads & Span-8 Pod.

•

- Color-coded racks
- Easy visual checks

Product Number: Beckman Coulter B85881; Thermo Fisher NC1392862

Platform: Biomek i-Series

Max Capacity Volume Range: 90 μ L

Pod and Tool Type: 300 μ L MC-96 Head, 1200 μ L MC-96 Head, Span-8 Pod

Tip Type: 90 μ L

Instrument Compatibility: Biomek i-Series

Package Quantity: 10 racks (96 per rack)

Autoclavable: No

Sterile: No

Color: Clear

Material(s): Polypropylene



Note: This product was over-purchased from Beckman Coulter™.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

