

BECKMAN COULTER BIOMEK I-SERIES TIPS 230UL NON-STERILE

Biomek® Pipette Tips from Beckman Coulter® are the only tips on the market that are designed, validated, and approved for use on Biomek® Liquid Handlers.

- Color-coded racks
- Easy visual checks

Tip definitions, pipetting techniques, and templates are built into Biomek® software. For best in total system performance, use only Beckman Coulter® validated pipette tips on Biomek® workstations.

Product Number: Beckman Coulter B85903; Thermo Fisher NC1392864

Platform: Biomek i-Series

Max Capacity Volume Range: 230 µL

Pod and Tool Type: 300µL MC-96 Head, 1200µL MC-96 Head, Span-8 Pod

Tip Type: 230 µL

Instrument Compatibility: Biomek i-Series

Package Quantity: 10 racks (96 per rack)

Autoclavable: No

Sterile: No

Color: Clear

Material(s): Polypropylene



Note: This product was over-purchased from Beckman Coulter™.