

BD LUER LOCK SYRINGE, 50ML

BD Luer Lock Syringe, 50mL, 40/Box 309653

- Feature a clear barrel with bold scale markings
- Tapered plunger rod for ease of aspiration, positive plunger rod stop
- Packaging is clearly labeled: free of latex, pyrogen, BPA, DEHP or PVC

Cat. #: BD 309653; Thermo Fisher Scientific 13-689-8

Qty./Pk: 40 / Pk

Brand: DB®

Disposable: Single-use


Graduations: 1 mL

Volume (Metric): 50 mL

Connection Type: Luer Lock

Sterility: Sterile

Needle Point Number: No Needle Included



Note: This product was over-stocked. It is 50% off Corning list price.