

BD SPINAL NEEDLES

Description:

The BD 405174 Spinal Needle is a disposable, 18-gauge x 3-inch sterile needle featuring a Quincke-style bevel and a color-coded pink hub for fast identification. Designed by BD (Becton, Dickinson and Company) for high-precision spinal procedures, this needle offers excellent tactile feedback, reliable CSF flow, and reduced patient trauma.

Material	Stainless Steel
Color	Pink
Needle Gauge	18 G
Includes	Key/slot arrangement of fitted stylet, cannula hubs, translucent window hub
Latex Free	Yes
Length (English)	3 in.
Product Type	Spinal Needle
Disposable	Single-use
Quantity	25/Pk.
Needle Point Style	Bevel
Sterility	Sterile