

## BIOTINYLATED ANTI-CD28 (DM64) RABBIT MAB

**Cat.#:** 28372

**Product Name:** Biotinylated Anti-CD28(DM64) Rabbit Monoclonal Antibody

**Synonyms:** CD28; Tp44

**Description:** Biotinylated Anti-CD28 antibody(DM64) Rabbit Monoclonal Antibody

**Background:** The protein encoded by this gene is essential for T-cell proliferation and survival; cytokine production; and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Applications:** ELISA; Flow Cyt

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Host Species:** Rabbit

**Isotype:** Rabbit IgG

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human CD28

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).