

BIOTINYLATED ANTI-GPRC5D (DM61) RABBIT MAB

Cat.#: 28366

Product Name: Biotinylated Anti-GPRC5D(DM61) Rabbit Monoclonal Antibody

Synonyms: GPRC5D

Description: Biotinylated Anti-GPRC5D antibody(DM61) Rabbit Monoclonal Antibody

Background: The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.

Applications: ELISA; Flow Cyt

Recommended Dilutions: IHC 1:200, ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Rabbit

Isotype: Rabbit IgG

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human GPRC5D

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).