

## BIOTINYLATED ANTI-GPRC5D (DM61) RABBIT MAB

**Cat.#:** 28366

**Product Name:** Biotinylated Anti-GPRC5D(DM61) Rabbit Monoclonal Antibody

**Synonyms:** GPRC5D

**Description:** Biotinylated Anti-GPRC5D antibody(DM61) Rabbit Monoclonal Antibody

**Background:** The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.

**Applications:** ELISA; Flow Cyt

**Recommended Dilutions:** IHC 1:200, ELISA 1:5000-10000; Flow Cyt 1:100

**Host Species:** Rabbit

**Isotype:** Rabbit IgG

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human GPRC5D

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).