

BIOTINYLATED ANTI-TROP2 (DM74) RABBIT MAB

Cat.#: 28390

Product Name: Biotinylated Anti-Trop2(DM74) Rabbit Monoclonal Antibody

Synonyms: TACSTD2; GA733-1; MIS1; TROP2

Description: Biotinylated Anti-Trop2 antibody(DM74) Rabbit Monoclonal Antibody

Background: This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Rabbit

Isotype: Rabbit IgG

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human Trop2

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).