

## BIOTINYLATED ANTI-TROP2 (DM75) RABBIT MAB

**Cat.#:** 28392

**Product Name:** Biotinylated Anti-Trop2(DM75) Rabbit Monoclonal Antibody

**Synonyms:** TACSTD2; GA733-1; MIS1; TROP2

**Description:** Biotinylated Anti-Trop2 antibody(DM75) Rabbit Monoclonal Antibody

**Background:** This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

**Applications:** ELISA; Flow Cyt

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Host Species:** Rabbit

**Isotype:** Rabbit IgG

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human Trop2

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).