

CD114 (DM172) RABBIT MAB

Cat.#: 28515

Product Name: Anti-CD114(DM172) Rabbit Monoclonal Antibody

Synonyms: CSF3R;CD114;GCSFR

Description: Anti-CD114 antibody(DM172) Rabbit Monoclonal Antibody

Background: The protein encoded by this gene is the receptor for colony stimulating factor 3; a cytokine that controls the production; differentiation; and function of granulocytes. The encoded protein; which is a member of the family of cytokine receptors; may also function in some cell surface adhesion or recognition processes. Alternatively spliced transcript variants have been described. Mutations in this gene are a cause of Kostmann syndrome; also known as severe congenital neutropenia.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Rabbit

Isotype: Rabbit IgG

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CD114

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-CD114 (DM172) on Expi293 cells transfected with human CD114 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).