

CORNING® 250 ML POLYCARBONATE ERLLENMEYER FLASK WITH VENT CAP

- Made from optically clear polycarbonate
- Ideal for shaker culture applications
- Sterilized by gamma radiation and nonpyrogenic
- RNase- and DNase-free

Cat. #: 431144

Qty./Pk: 1 / Pk

Qty./Cs: 50 / Cs

Brand: Corning®

Neck Style: Wide Mouth

Cap Style: Vented

Cap Material: Polypropylene

Cap Color: Orange

Neck O.D.: 31 mm (approx.)

Sterile: Yes

Feature: Disposable

Capacity: 250 mL

Flask Style: Erlenmeyer

Top Style: Threaded

Baffles: No

Graduation Interval: 25 mL



Note: This product was over-purchased from Corning. It is 70% off Corning list price.