

## **CORNING® 500 ML PP CENTRIFUGE TUBES WITH PLUG SEAL CAP, STERILE, 6/PACK, 36/CASE**

Corning® 500 mL polypropylene (PP) centrifuge tubes are ideal for applications requiring large volume centrifugation. Please check our Technical Information Listing for information on the physical and chemical properties of Corning plastic centrifuge tubes as well as important suggestions for safe use.

- Support cushions (Catalog No. 431124 for 500) must be used with this product unless the centrifuge rotor has appropriate V-bottom cushions.
- High Density Polyethylene plug seal caps
- Maximum RCF of 6,000xg.
- Tubes are sterile and nonpyrogenic

Product Number: 431123

Qty./Pk: 6 / Pk

Qty./Cs: 36 / Cs

Brand: Corning®

Capacity: 500 mL

Tube Style: Centrifuge

Top Style: Threaded

Maximum RCF: 6000 xg

Graduations: Yes

Color: Natural

Cap Style: Plug seal

Cap Material: High density polyethylene

Cap Color: Orange

Sterile: Yes



**Note:** This product was over-purchased from Corning. It is 70% off Corning list price.